

VIBE BATTERY CARE

The **VIBE** contains a **sealed lithium battery**, designed for safety and durability, which isn't user-replaceable. However, with good charging habits, most users experience **well beyond 5 years** of reliable use. Here's how to maximize its life span:



Best Practices to Extend Battery Life

1 Avoid overcharging.

Once the battery indicator shows full (solid red light), unplug it. Continuous overcharging shortens lithium cell life.

2 Don't let it drain completely.

Try to recharge when the battery level drops below 20%, rather than waiting for it to go flat.

3 Use only standard 5V chargers.

Avoid fast chargers or car USB ports, which can push higher voltages and strain the battery.

4 Store at moderate temperature.

Keep your device in a cool, dry environment. Avoid leaving it in direct sunlight, hot cars, or cold areas.

5 Regular use helps.

Running the device regularly (a few times per week) keeps the battery cells active and stable.



Battery & Technology Notes

- **Runtime:** 4–6 hours on a full charge.
- **Recharge time:** 1-2 hours with a standard 5V adapter.
- The sealed battery design is what allows the VIBE to remain compact, lightweight, and free of vibration or noise.
- Resona is always developing new wellness technology if a **replaceable-battery version** becomes available in the future, current owners will be notified via our newsletter

Many of our early users have been running their original units for **over four years** with no noticeable drop in charge time, simply by following the care steps above.

FDA Disclaimer: This product is intended for general wellness purposes only and is not intended to diagnose, treat, cure, or prevent any disease.